

PRODUCT SPECIFICATION

Product PN FI148AV29N200A00

Mod. 984 c

Description Tubular Filter with PP frame and NY mesh 200 um

Rev. 03

Tubular filter PP frame + 200 um PA mesh



PRODUCT DESCRIPTION	Non-sterile tubular filter with mesh Ny 200µm Product provided in bulk packs for further manufacturing
INTENDED USE / APPLICATION	Transfusion
MATERIALS	Mesh: Polyamide 200 um Housing: Polypropylene Regulatory Documentation Required: <ul style="list-style-type: none"> <input checked="" type="checkbox"/> Biocompatibility according USP VI <input checked="" type="checkbox"/> Rohs, directive 2002/95/CE <input checked="" type="checkbox"/> BSE/TSE, directive 2003/32/CE <input checked="" type="checkbox"/> DEHP plasticizer Free and latex free <input checked="" type="checkbox"/> Reach 1907/2006/CE (hazardous substances regulation) <input checked="" type="checkbox"/> Conflict minerals
PRODUCT SHELF LIFE	5 years
STERILIZATION	STERILE : <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No SUITABLE FOR STERILIZATION : <input checked="" type="checkbox"/> Ethylene oxide
PACKAGING AND LABELLING	8.000 pcs packed in two PE bags. Bags are separately closed. The first bar-code label is between the 2 bags. The second bar-code label is stuck outside the box. Each bag is labeled with the following traceability information: <ul style="list-style-type: none"> - Quantity - Product description - Product date - Lot number (OL and batch number to trace back to raw materials used) - Operator code Different lot of goods in one shipment are packed in a manner to prevent mix-ups. Different lots in one box are separately closed and separately labeled to prevent mix-ups.

PRODUCT SPECIFICATION

Product PN	FI148AV29N200A00	Mod. 984 c
Description	Tubular Filter with PP frame and NY mesh 200 um	Rev. 03

CERTIFICATE OF COMPLIANCE	Conformity declaration is printed on every invoice and Certificate is according to UNI EN 10204 type 2.1 The Quality management system in compliance with ISO 9001, ISO 13485.													
ACCEPTABLE QUALITY LEVEL	AQL (see table below)													
VISUAL REQUIREMENTS	<p><i>Visual acceptance requirements apply when inspected under below conditions:</i></p> <p>Instrument inspection with naked eye and light source Distance of 300-450 mm. Timings: 5 sec per unit</p> <table> <tr> <th colspan="2">Acceptance Requirement</th><th>AQL</th><th>Sampling Plan</th></tr> <tr> <td>1</td><td> 1. Short molding 2. Mesh threads protruding from plastic ring (filter top area) 3. Short mesh, detached from the plastic support. 4. Broken mesh 5. Open mesh out of specification (uneven mesh). 6. Missing mesh </td><td>0.4</td><td rowspan="2">ISO 2859 part. 1 1st Level</td></tr> <tr> <td>2</td><td> 7. Mesh protruding out of the rib 8. Mesh welding out of the rib 9. Flashes > 0.3 mm. 10. Gate scar > 0.3 mm. 11. Loose particulate matter > 0.1 mm² 12. Plastic infiltration on filter mesh > 2nd mesh opening 13. Dents/Scratches with plastic residual 14. Shrinkage on plastic ring (filter top area) 15. Burning marks 16. Dirty of oil, grease 17. Embedded particulate matter > 0.2 mm² </td><td>1.5</td></tr> </table> <p>*Embedded Particulate Matter, loose contamination and air bubbles: according to Dirt Estimation Chart (Tappi Standard). Characteristics listed above are statistically inspected during the manufacturing process.</p>			Acceptance Requirement		AQL	Sampling Plan	1	1. Short molding 2. Mesh threads protruding from plastic ring (filter top area) 3. Short mesh, detached from the plastic support. 4. Broken mesh 5. Open mesh out of specification (uneven mesh). 6. Missing mesh	0.4	ISO 2859 part. 1 1 st Level	2	7. Mesh protruding out of the rib 8. Mesh welding out of the rib 9. Flashes > 0.3 mm. 10. Gate scar > 0.3 mm. 11. Loose particulate matter > 0.1 mm ² 12. Plastic infiltration on filter mesh > 2nd mesh opening 13. Dents/Scratches with plastic residual 14. Shrinkage on plastic ring (filter top area) 15. Burning marks 16. Dirty of oil, grease 17. Embedded particulate matter > 0.2 mm ²	1.5
Acceptance Requirement		AQL	Sampling Plan											
1	1. Short molding 2. Mesh threads protruding from plastic ring (filter top area) 3. Short mesh, detached from the plastic support. 4. Broken mesh 5. Open mesh out of specification (uneven mesh). 6. Missing mesh	0.4	ISO 2859 part. 1 1 st Level											
2	7. Mesh protruding out of the rib 8. Mesh welding out of the rib 9. Flashes > 0.3 mm. 10. Gate scar > 0.3 mm. 11. Loose particulate matter > 0.1 mm ² 12. Plastic infiltration on filter mesh > 2nd mesh opening 13. Dents/Scratches with plastic residual 14. Shrinkage on plastic ring (filter top area) 15. Burning marks 16. Dirty of oil, grease 17. Embedded particulate matter > 0.2 mm ²	1.5												

PRODUCT SPECIFICATION

Product PN FI148AV29N200A00

Mod. 984 c

Description Tubular Filter with PP frame and NY mesh 200 um

Rev. 03

DIMENSIONAL REQUIREMENTS

Acceptance Requirement		AQL	Sampling Plan
Total height	40.85 ± 0.15 mm	0,1	ISO 2859 part. 1 S3 Special Inspection Levels
External diameter	15 ± 0.15 mm		
Base height	4 ± 0.10 mm		

This material specification describes the properties of product above indicated.

This document contains general requirements, material description, drawing references, defect specification, biological material requirements.

REVISIONS AND APPROVALS:

DATE	REV.	REASON FOR CHANGE	ISSUED AND CONTROLLED BY: (name /function and signature)	APPROVED BY: (name /function and signature)
20/05/2019	00	Validation of GVS Microfiltration	Paolo Carati (RP) <i>Paolo Carati</i>	Sara Zaccheroni (AQ) <i>Sara Zaccheroni</i>

Customer Approval:

We accept this material specification as a part of the agreed terms of delivery

Company name _____

Approved by: _____
(Name, Function) (Signature)

Date _____
(Company stamp)

Please send back this document signed for approval. If we will not receive this specification signed , we consider the first order placed as implicit approval.