




Troubleshooting




Problem	Cause	Solution
High background 	Insufficient blocking	- Block the membrane for at least one hour with skim milk or BSA depending on the type of protein.
	Blocking is not compatible	- For phosphorylated protein detection, milk should not be used (milk and casein are phospho-protein rich).
	High concentration of primary antibody	- Reduce and/or optimize the concentration of the primary antibody.
	High concentration of secondary antibody	- Reduce and/or optimize the concentration of the secondary antibody.
	Long exposure time with chemiluminescent substrate	- Decrease the membrane incubation time in the substrate. - Decrease the blot detection time, e.g. 120 or 60 seconds instead of 180 seconds.
	Too sensitive Detection reagents	- In the case of using a strong chemiluminescent substrate for the detection step, in order to reduce the background while having the lanes, dilute the substrate with PBS (or other buffers you are using for washing without detergent) with a 2: 1 proportion. - Decrease the quantity of the substrate (put less volume of each solution but still a 1:1 proportion). Follow the substrate supplier's instructions.
	The membrane has dried out	- Never let the membrane dry out during the Western Blotting workflow in the case of hydrophobic membranes.
	Membrane with an excess of chemiluminescence substrates	- Before imaging, carefully drip out the excess substrate from the membrane.
	The membrane was in contact with the container in the equilibrating or washing step	- Keep the protein side of the blot upward and in contact with the solution. - Wash the membrane at least three times, each for a minimum of 5 minutes.
	Insufficient washing	- Insufficient volume of washing buffer. - Insufficient concentration of the detergent in the buffer. Use PBS containing >0.05% Tween 20.

Problem	Cause	Solution
<p data-bbox="113 208 568 264">Patchy or Uneven black Spots on the Blot</p>  	<p data-bbox="568 253 1023 275">The membrane got dried during the detection stage</p> <hr/> <p data-bbox="568 383 1023 405">Contaminated membrane</p> <hr/> <p data-bbox="568 456 1023 524">The secondary antibody is aggregated µm RC, GVS).</p> <hr/> <p data-bbox="568 584 1023 607">Antibodies are binding to the blocking buffer</p> <hr/> <p data-bbox="568 669 1023 692">Insufficient blocking</p> <hr/> <p data-bbox="568 757 1023 779">Contaminated equipment</p> <hr/> <p data-bbox="568 862 1023 884">Contaminated reagents</p> <hr/> <p data-bbox="568 927 1023 949">No more activity of reagents</p>	<p data-bbox="1023 181 1473 226">- Incubate the membrane in the detection substrate according to the producer's recommendation.</p> <p data-bbox="1023 237 1473 349">- To avoid drying during the detection step, place the blot in between the development folder and gently smooth out the air bubbles. Be careful not to dry the membrane.</p> <hr/> <p data-bbox="1023 356 1473 378">- Use clean buffers in each step.</p> <p data-bbox="1023 385 1473 430">- Be careful about the residues of the gel on the membrane.</p> <hr/> <p data-bbox="1023 441 1473 497">- Filter the secondary antibody using a 0.22 µm low protein binding filter (e.g. SEPARA 0.2</p> <hr/> <p data-bbox="1023 539 1473 618">- Filter the blocking solution antibody using a 0.22 µm filter to remove undissolved blocking agent or any contamination.</p> <p data-bbox="1023 629 1473 651">- Use an alternative blocking agent.</p> <p data-bbox="1023 658 1473 703">- Block the membrane for at least one hour with skim milk or BSA depending on the type of protein.</p> <p data-bbox="1023 714 1473 770">- Do not use metallic forceps, which can cause speckling due to rust contamination.</p> <p data-bbox="1023 781 1473 826">- Use clean plastic containers to avoid any type of cross-reaction.</p> <p data-bbox="1023 837 1473 882">Make sure that all solutions, buffers, and chemiluminescence substrates are clean</p> <p data-bbox="1023 893 1473 916">- Avoid repetitive freezing/thawing of the proteins.</p> <hr/> <p data-bbox="1023 927 1473 949">- Check the datasheet and expiration date of reagents.</p>
<p data-bbox="113 947 568 1003">No band (including protein marker) or fade bands</p> 	<p data-bbox="568 1003 1023 1025">Too strong washing</p> <hr/> <p data-bbox="568 1205 1023 1227">No transfer</p>	<p data-bbox="1023 969 1473 992">- Try to reduce repetition and/or time of washing.</p> <p data-bbox="1023 1003 1473 1059">- Reduce the concentration of detergent (e.g. Tween 20) in the washing buffer.</p> <hr/> <p data-bbox="1023 1137 1473 1182">- Check the transfer efficacy using Ponceau S as a reversible stain on the membrane.</p> <p data-bbox="1023 1193 1473 1238">- Check the transfer efficacy using coomassie Blue on the gel.</p> <p data-bbox="1023 1249 1473 1305">- Check if the "sandwich" (sponge, filter paper, gel, and membrane) is assembled and placed correctly.</p>
<p data-bbox="113 1395 568 1451">No band (but protein marker present)</p> 	<p data-bbox="568 1429 1023 1451">Wrong antibodies</p> <hr/> <p data-bbox="568 1653 1023 1675">The activity of antibodies is reduced</p>	<p data-bbox="1023 1384 1473 1429">- Check the compatibility of the primary antibody with your target protein.</p> <p data-bbox="1023 1440 1473 1496">- Check the compatibility of the primary antibody with the secondary antibody.</p> <hr/> <p data-bbox="1023 1630 1473 1653">- Try a fresh antibody.</p> <p data-bbox="1023 1664 1473 1686">- Check the storage condition.</p> <p data-bbox="1023 1697 1473 1720">- Avoid repetitive freezing/thawing of the antibodies.</p>
<p data-bbox="113 1843 568 1888">Uneven white spots</p> 	<p data-bbox="568 1955 1023 1977">Presence of air bubbles during the transfer</p>	<p data-bbox="1023 1933 1473 2011">- Air bubbles between gel and membrane can cause white spots. Use a proper roller to smooth out all the bubbles before starting the transfer.</p>

Problem	Cause	Solution
---------	-------	----------

Smear bands



Too much sample is loaded

- Dilute more the loading sample.

Poor quality of the sample

- Contaminated, degraded, or denatured sample.

Poor gel preparation (uneven gel)

- Fill the cassette to the volume before the gel starts polymerization.


White/hollow bands



High concentration of reagents

- Decrease the concentration of primary/secondary antibodies or use less protein.

Bands smile within the gel lanes



The running condition was too fast

- Reduce voltage during electrophoresis.

Migration was too hot

- Run the gel in the cold room.

Broad or misshapen bands



Gel electrophoresis problems

- Poor gel polymerization.
- Inappropriate running condition.