HI-FLO IV Sets
HI-FLO Extension Sets
Speedflow® IV Filters
EPI-Baby IV Filters

Instructions for Use
Istruzioni d'uso
Gebrauchsanweisung
Mode d'emploi
Instrucciones de uso

 **GVS S.p.A**
Via Roma 48/50
40069 Zola Predosa (BO) - Italy
Tel.: +39 051 6176311
Fax: +39 051 6176200
E-mail: event-notification@gvs.it
Web: www.gvs.com

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Warning: Non-compliance with these instructions for use might lead to user/patient injuries or other unforeseen side effect.

Warning: Any serious incident that has occurred in relation to the device must be reported to GVS S.p.A. and the competent authority of the Member State in which the user and/or patient is established.

1. ⚠ Warnings

- Warning:** Do not use if packaging is open or appears damaged.
Warning: Do not use if product appears damaged.
Warning: Do not use if protective caps, if present, are disconnected or lost.
Warning: Do not use after the indicated expiration date.
Warning: Do not re-sterilize.
Warning: Disposable and single patient use. Do not re-use. Reuse involves the risk of infection to the patient and/or user and may compromise the functionality of the device.
Warning: If the extension is intended to be connected with a terminal filter, it must not be used with infusions or infused with liquids which may clog the filter.
Warning: Use the device only when compatibility with the drugs/solutions is known. It's recommended to test or verify the compatibility between drugs/substance and device prior to its use.
Warning: Non-compliance with these IFUs may lead to users/patients injuries or other unforeseeable risk.
Warning: The manufacturer disclaims any liability for patient injury caused by improper use of the device and / or other than that provided in the instructions for use/label or resulting from use by unqualified or untrained.
Warning: to ensure the correct storage and handling keep the devices in dry conditions and do not expose the device to direct sunlight. Do not use the device if you suspect that it has been stored in wet conditions or under direct sunlight.
Warning: use the device according to clinical aseptic procedures and hospital protocol.
Warning: do not use isopropyl alcohol to clean the device and do not connect the device to parts cleaned with such detergent. Only devices with REF. GMD019A00S are isopropyl alcohol resistant.
Warning: it's recommended to wear gloves when priming the infusion set, especially when dangerous drugs/substances are handled.
Warning: Do not use if sensitivity between patient and device materials are known;
Warning: Do not use if incompatibility with fluids being filtered is known;
Warning: If the Y-site is after the filter, do not infuse fluids through the Y-site if filtration of these fluids/liquid mixture is required;
Warning: Do not use if incompatibility occurs between the infused drugs and the device materials.

2. ⚠ Precautions and Recommendations

- Precaution:** If flow rate decreases, replace IV Set.
Precaution: The Infusion Controllers are used for administration in gravity, while the Infusion extensions and filters may be used in infusion devices that either work by gravity or connected to a peristaltic pump or syringe. Filters may also be used as syringe filters or connected to other infusion sets. The intended application as per gravity or with a pump is indicated on the label.
Precaution: For Extensions with vented filters and vented filters of 0,2, 0,2+ or 1,2 µm pore sizes do not exceed the pressure of 3,2 bar (as indicated on the label). For Extensions with non-vented filters and non-vented filters of 0,2 µm pore sizes do not exceed the pressure of 8,0 bar (as indicated on the label).
Precaution: Devices/Filters with 1,2 µm pore sizes shall not be used for more than 24 hours. Devices /Filters with 0,2 µm pore sizes shall not be used for more than 96 hours. Devices /filters with 0,2+ µm pore sizes shall not be used for more than 120 hours.
Precaution: Do not excessively pull the tubes to prevent damages or disconnections.
Precaution: Do not use flow regulators near the OPEN position because the flow would be too fast (more than 2,5 lt/h).
Recommendation: When the flow rate is lower than 30 ml/h, the recommendation is to count the number of drops/minute and to adjust the flow by raising or lowering the bottle or by setting the flow regulator/roller to the required flow rate.
Recommendation: Easydrop®/Eurodrop® is not a measurement instrument; flow rate is produced by gravity. Many factors, such as the patient position, his venous pressure and the bottle height can affect the flow rate.
Recommendation: Adherence with local protocols regarding aseptic technique is recommended.
Recommendation: The whole system and connections has to be perfectly locked before use.
Recommendation: Patients must be informed about any contraindication, residual risk, side-effect, warning, precaution and measure.

3. Device Description

Device for infusion, with different configurations: Infusion Controllers, Infusion extensions, Infusion filters.
 The products are delivered in sterile conditions and reach their intended purpose by gravity or activated by peristaltic pumps, depending on the type of device (refer to par. 4).

Can include a filter with hydrophilic PES membrane 0,2/0,2+/1,2 or 5,0 µm pore size for particles retention (printed in black on the housing as "0,2 micron", "0,2 positive", "1,2 micron" or "5,0 micron" respectively). Filter may be vented or non-vented. The vented versions have a PTFE membrane 0,03 µm pore size.
 HI-FLO IV Sets and HI-FLO Extension Sets may include an Easydrop®/Eurodrop® flow regulator single or double scale, to regulate the flow rate. The single scale version has a general purpose scale that can be used for any kind of solution. The double scale version has one scale for liquids having density lower than 10% and a second scale for liquids with higher density.
 HI-FLO IV Sets and HI-FLO Extension Sets may also include:

- Vented IV Drip Chamber (only on HI-FLO IV sets)
- Roller regulator
- Y-Injection Point
- Clamp

4. Intended Use

Devices are intended for infusion. They can be used for intravascular or neuraxial applications. It's up to the user to choose the applicable infusion set on the basis of the specific application, the liquid to be infused and the established clinical practice.

- Infusion Controllers (HI-FLO IV Sets) are for gravity feed infusion only;
- Infusion Extensions (HI-FLO Extension Sets) with EasyDrop®/EuroDrop® flow regulator are infusion extension lines for use by gravity;
- Infusion Extensions (HI-FLO Extension Sets) without EasyDrop®/EuroDrop® flow regulator are infusion extension lines that can be used in combination with an infusion set or an infusion device (pump or syringe);
- Infusion filters (Speedflow® IV Filters, EPI-Baby IV Filters) can be used in combination with an infusion set or an infusion device (pump or syringe).

Sets equipped with IV filters (Adult, Pediatric or Neonatal) are used specifically to eliminate air bubbles and to retain particles, bacteria (0,2 and 0,2+ µm filters), endotoxins (0,2+ µm filter), fungi and microorganisms (1,2 µm filters), particles (5,0 µm filters).
 Infusion Devices equipped with a flow regulator allow to set approximately the quantities infused.
 In case of devices intended for use with pressure infusion equipment the maximum working pressure is indicated on the label.
 Infusion Devices are single-use devices that can be used applications that last:

- Devices with 0.2 µm filter: up to 96 hours;
- Devices with 0,2+ µm filter: up to 120 hours;
- Devices with 1.2 µm filter: up to 24 hours.
- Devices without filter don't have a specific duration but are single use and have to be disposed after each therapy.

Infusion Devices are supplied sterile.
 Infusion Devices may be used for administration of the following substances in accordance with their manufacturers' instructions:

- Drugs
- Anesthetics
- Nutritional Supplements
- Saline solutions

Based on the tests carried out, GVS S.p.A. can conclude that there is no significant absorption by Infusion Devices of the following drugs/solutions, no further compatibility tests have been performed:

- Lidocaine HCL (Anesthetic)
- Bupivacaine (Anesthetic)
- Piperacillin Sodium (Antimicrobial)
- Paclitaxel (Chemotherapeutical)
- Fluorouracil (Chemotherapeutical)
- Etoposide (Chemotherapeutical)
- Dobutamine (Cardiovascular)
- Sodium Citrate (Anticoagulant)
- Insulin (Protein Fraction)
- Mannitol (Diuretic)
- Common Saline Solution

The infusion systems should only be supervised and used by qualified healthcare professionals. Use environment is the medical one, typically hospital, nursing homes and homes under supervision from a healthcare professional.
 Products are designed to be used with adults, pediatric and neonatal patients as specified on each device. Weight, Gender, Nationality are not relevant.

4.1 ⚠ Contraindications

- Configurations of sets/filters with 0,2/0,2+ µm pore size cannot be used for infusion of hyper alimentation liquids, lipids, blood and/or blood derivatives;
- Configurations of sets/filters with 1,2 µm pore size cannot be used for infusion of blood and/or blood derivatives;
- Do not use the sets/filters with 5,0 µm pore size in infusion therapy without further bacterial retention filter.

5. Delivery Status and Condition of Use

Devices are delivered sterile, sterilized using ethylene oxide.
 Date of manufacture and use-by date are indicated on the label.

6. Device Combination

Infusion Devices can be connected to all administration devices and components such as needles and syringes commercially available (compliant with ISO 80369-7).

Extensions and Filters for use with pressure devices can be connected to infusion pumps using an ISO 80369-7 connector, given that the maximum pressure of the pump is equal or lower than the maximum pressure indicated on the device label. Depending on the type of pump (e.g. elastomeric pump vs peristaltic pump) and on the presence of an ISO 80369-7 connection, devices may be connected directly to the pump or not. In the last case the user has to use connection devices (e.g. tubes) compliant with ISO 80369-7.

7. Instructions for use

7.1 Infusion Controllers (HI-FLO IV Sets)

1. After removing the protective cap, connect the device to the bottle or to the infusion container, ensuring that the clamp is in the closed position and the flow regulator (or the roller) off.
2. Hang the bottle, placing the drip chamber about 80 cm above the patient's middle axillary line.
3. To eliminate the air from the set, open the clamp and set the flow regulator (or the roller) to the OPEN position. If present, place the filter in a vertical position with the arrow pointing up. Wait for the complete evacuation of air. Then close the clamp.
4. Check that the system is connected and tighten correctly.
5. After aseptically removing the protective cover from the outlet connector, connect the set to the venous puncture device.
6. Set the flow regulator, if present, to the desired flow rate.
7. Check the flow rate by counting the number of drops/minute. Raise or lower the infusion container to achieve the required flow rate.

Precaution: The following table has been set out for NaCl 0,9% solution. If running an IV solution of a higher viscosity than NaCl 0,9%, the flow rate must be adjusted by counting drops. For solutions as e.g. glucose 40% the second scale can be used for orientation.

Infusion Time [h]	100 [ml]		250 [ml]		500 [ml]		1000 [ml]	
	Flow Rate [ml/h]	Drops /min	Flow Rate [ml/h]	Drops /min	Flow Rate [ml/h]	Drops /min	Flow Rate [ml/h]	Drops /min
1	100	33	250	83	500	167	1000	333
2	50	17	125	42	250	83	500	167
3	33	11	83	28	167	56	333	111
5	20	7	50	17	100	33	200	67
7	14	5	36	12	71	24	143	48
9	11	4	28	9	56	19	111	37
12	-	-	21	7	42	14	83	28
18	-	-	14	5	28	9	56	19
24	-	-	10	3	21	7	42	14

7.2 Infusion Extensions (HI-FLO Extensions Set) with Easydrop®/Eurodrop®

1. Remove the covers from the inlet and outlet connections and aseptically connect the set to the IV administration line.
2. Connect the IV administration device to the infusion container, making sure the flow regulator is in the OFF position.
3. Hang the bottle so that the drip chamber is approximately 80 cm above the patient's middle axillary line.
4. To eliminate the air from the system, adjust the flow regulator to the OPEN position, aseptically remove the protective cap of the male luer lock to let the liquid come out from the lower end, then return the flow regulator to the OFF position.
5. Check that the system is connected and tighten correctly.
6. Connect the outlet connector of the infusion line to the venous puncture device.
7. Set the flow regulator to the desired flow rate.
8. Check the flow rate by counting the number of drops/minute. Raise or lower the bottle to achieve the required flow rate.

7.3 Infusion Extensions (HI-FLO Extension Sets) with filter

1. Remove the covers from the inlet and outlet connections and aseptically connect the set to the IV administration line.
2. Ventilate the line by keeping the filter in a vertical position with the arrow pointing up to prevent air bubbles from partially blocking the surface of the filter. The filter must be filled from bottom to top. Upon completion of the priming the filter will work in all positions.
3. Connect the line to the intended point.
4. Check that the system is connected and tighten correctly.

7.4 Infusion Filters (Speedflow® IV Filters, EPI-Baby IV Filters)

1. Remove the covers from the inlet and outlet connections and aseptically connect the filter to the IV administration line.
2. Proceed with the priming. Place the filter in vertical position with direction arrow upward. Filter must be filled from bottom to top to avoid air bubbles blocking partially the filter surface.
3. Check that the system is connected and tightened correctly.

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