



PRODUCT SPECIFICATION

CE 0051

Product PN	GMD001A30S GMD002A30S	Mod. 984e-ext
Description	Speedflow® IV filter Adult 0.2/1.2 µm Non-Vented	Rev. 07

Speedflow® filter Adult 0.2/1.2 µm Non-Vented



PRODUCT DESCRIPTION	<p>Non-vented filter for and anaesthesia.</p> <p>The device is sterile and can be used in combination with an infusion set or an infusion device (pump or syringe). The connections are a female luer lock with vented cap inlet and a rotating male luer lock with vented cap outlet, compliant with ISO 80369-7.</p> <p>The device is provided with a non-vented Speedflow® Adult filter, with an hydrophilic PES membrane with 0.2/1.2 µm pore size for particles retention.</p> <p>The filter size (LxWxH) is 30x67x7,2 mm.</p>	
CONFIGURATIONS	GMD001A30S	0.2 µm
	GMD002A30S	1.2 µm
MANUFACTURER NAME	<p>GVS SpA Via Roma, 50 40069 Zola Predosa (BO) – ITALY Phone: +39.051.6176311 – Fax: +39.051.6176200 e-mail: gvs@gvs.it – website: www.gvs.com</p>	
INTENDED USE / APPLICATION	<p>Infusion of Drugs (including chemotherapy), Anesthetics, Nutritional Solutions and Saline Solutions.</p> <p>The filter can be connected via lock connectors to all commercially available infusion administration sets compliant with ISO 80369-7.</p> <p>The filter can be used in combination with an infusion set or an infusion device (pump or syringe).</p> <p>The device with 0.2 µm filter is intended specifically to eliminate air bubbles and to retain particles and bacteria, while the device with 1.2 µm filter is intended specifically to eliminate air bubbles and to retain fungi and microorganisms.</p> <p>The device is a single-use device that can be used applications that last:</p> <ul style="list-style-type: none">▪ Devices with 0.2 µm filter: up to 96 hours;▪ Devices with 1.2 µm filter: up to 24 hours. <p>The device has to be disposed after each therapy.</p> <p>Devices with 0.2 µm filter cannot be used for infusion of hyper alimententation liquids, lipids, blood and/or blood derivatives.</p> <p>Devices with 1.2 µm filter cannot be used for infusion of blood and/or blood derivatives.</p> <p>The device should only be supervised and used by qualified healthcare personnel.</p>	
CLASS OF THE PRODUCT	<p>Class IIa – sterile – single use</p> <p>Rule 3 Annex VIII Regulation (UE) 2017/745</p>	
EMDN	A04010101	WITHDRAWAL AND ADMINISTRATION FILTERS (0.2-1.2 MICRONS)
MATERIALS	<p>Filter media: PES</p> <p>Filter Housing: PP, MABS</p> <p>Caps: HDPE</p> <p>Regulatory Compliance</p> <ul style="list-style-type: none">▪ Biocompatibility according to ISO 10993-1▪ Rohs directive 2011/65/UE▪ DEHP plasticizer Free▪ Latex free▪ Reach 1907/2006/CE (hazardous substances regulation)▪ Reg. 1272/2008/CE (classification, labelling and packaging of substances and mixtures)	



PRODUCT SPECIFICATION



Product PN	GMD001A30S GMD002A30S	Mod. 984e-ext Rev. 07
Description	Speedflow® IV filter Adult 0.2/1.2 µm Non-Vented	


PRODUCT CHARACTERISTIC	Physical/Mechanical: Dimensions (LxWxH): 30x67x7,2 mm Total internal volume of the set: 1.2 ml Weight: 7.86 g Input/output connectors: Female Luer Lock with vented cap inlet and rotating male Luer Lock with vented cap outlet compliant with ISO 80369-6 Operating temperature Range: From 5 °C to 40 °C Storage temperature Range: From 0 °C to 40 °C Maximum applicable pressure: 8.0 bar Flow Rate: 0.2 µm: ≥ 15 ml/min @ 80 cm (31.5 in) water head pressure 1.2 µm: ≥ 90 ml/min @ 80 cm (31.5 in) water head pressure Bubble point: 0.2 µm: 3.7 ÷ 4,8 bar 1.2 µm: 0.7 ÷ 1,0 bar
	Biological: Biocompatibility: Compliant with ISO 10993-1 Features: Type of administration: gravity / pressure Duration of the application: 0.2 µm up to 96 hours 1.2 µm up to 24 hours Filter: Speedflow® Adult 0.2/1.2 µm non-vented Filter pore size: 0.2/1.2 µm Filter internal volume: 1.2 ml Filter media: hydrophilic PES membrane
PRODUCT SHELF LIFE	5 years.
STERILIZATION	Sterile: Yes –Ethylene Oxide (EtO)
	Suitable for Resterilization: No
APPLICABLE STANDARDS AND REGULATIONS	Product Certification: CE mark
	Applicable Standards and Technical Regulations: EN 556-1 Sterilization of medical devices - Requirements for medical devices to be designated "STERILE" - Part 1: Requirements for terminally sterilized medical devices EN ISO 8536-4 Infusion equipment for medical use — Part 4: Infusion sets for single use, gravity feed EN ISO 8536-9 Infusion equipment for medical use — Part 9: Fluid lines for single use with pressure infusion equipment EN ISO 8536-11 Infusion equipment for medical use — Part 11: Infusion filters for single use with pressure infusion equipment EN ISO 10993-1 Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process EN ISO 10993-4 Biological evaluation of medical devices — Part 4: Selection of tests for interactions with blood EN ISO 10993-5 Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity EN ISO 10993-7 Biological evaluation of medical devices - Part 7: Ethylene oxide sterilization residuals EN ISO 10993-10 Biological evaluation of medical devices — Part 10: Tests for irritation and skin sensitization EN ISO 10993-11 Biological evaluation of medical devices -- Part 11: Tests for systemic toxicity EN ISO 10993-18 Biological evaluation of medical devices - Part 18: Chemical characterization of materials EN ISO 10993-23 Biological evaluation of medical devices — Part 23: Tests for irritation EN ISO 11135 Sterilization of health-care products — Ethylene oxide — Requirements for the development, validation and routine control of a sterilization process for medical devices EN ISO 11607-1 Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems EN ISO 11607-2 Packaging for terminally sterilized medical devices - Part 2: Validation requirements for forming, sealing and assembly processes

Product PN GMD001A30S GMD002A30S

Mod. 984e-ext

Description Speedflow® IV filter Adult 0.2/1.2 µm Non-Vented

Rev. 07

	<p>EN ISO 13485 Medical devices - Quality management systems - Requirements for regulatory purposes</p> <p>EN ISO 14971 Medical devices - Application of risk management to medical devices</p> <p>EN ISO 15223-1 Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements</p> <p>EN ISO 20417 Medical devices - Information to be supplied by the manufacturer</p> <p>IEC 62366-1 Medical devices – Part 1: Application of usability engineering to medical devices</p> <p>EN ISO 80369-1 Small-bore connectors for liquids and gases in healthcare applications — Part 1: General requirements</p> <p>EN ISO 80369-7 Small-bore connectors for liquids and gases in healthcare applications — Part 7: Connectors for intravascular or hypodermic applications</p> <p>EN ISO 80369-20 Small-bore connectors for liquids and gases in healthcare applications — Part 20: Common test methods</p>								
INSTRUCTIONS FOR USE	Available in: English / Italian / German / French / Spanish								
PACKAGING	<p>Primary Packaging: Devices are individually packed and label in medical paper pouches. Pouch Size: 100 X 145 mm</p> <p>Secondary Packaging: Pouches are placed inside a microperforated bag.</p> <p>Tertiary Packaging: Carton Box. Box Size: 60 x 40 x 20 cm Box Weight: 3 Kg Devices per box: 300</p>								
CERTIFICATIONS	<p>ISO 9001:2015</p> <p>ISO 13485:2016</p> <p>CE Certificate (2017/745/UE)</p>								
DRAWING	<div>  </div> <table border="1"> <thead> <tr> <th>ID</th><th>Description</th></tr> </thead> <tbody> <tr> <td>1</td><td>Vented male rotating cap for female Luer Lock</td></tr> <tr> <td>2</td><td>'1' model: Speedflow Adult 0.2 µm IV filter non-vented '2' model: Speedflow Adult 1.2 µm IV filter non-vented</td></tr> <tr> <td>3</td><td>Vented protecting cap for rotating male Luer Lock</td></tr> </tbody> </table>	ID	Description	1	Vented male rotating cap for female Luer Lock	2	'1' model: Speedflow Adult 0.2 µm IV filter non-vented '2' model: Speedflow Adult 1.2 µm IV filter non-vented	3	Vented protecting cap for rotating male Luer Lock
ID	Description								
1	Vented male rotating cap for female Luer Lock								
2	'1' model: Speedflow Adult 0.2 µm IV filter non-vented '2' model: Speedflow Adult 1.2 µm IV filter non-vented								
3	Vented protecting cap for rotating male Luer Lock								





PRODUCT SPECIFICATION

CE 0051

Product PN	GMD001A30S GMD002A30S	Mod. 984e-ext Rev. 07
Description	Speedflow® IV filter Adult 0.2/1.2 µm Non-Vented	

REVISIONS AND APPROVALS:

DATE	REV.	REASON FOR CHANGE	ISSUED BY: (name/function/signature)	VERIFIED BY: (name/function/signature)	APPROVED BY: (name/function/signature)
26/02/2023	00	First issue	Barbara Palmieri RA Specialist 	Barbara Finessi QA Manage 	Luca Zanini VP Healthcare and Lifesience 