

Product PN RS051- RS052- RS053- RS054- RS055- RS056- RS057- RS058 - RS062 – RS078 – RS079

Mod. 984 c

Description Baby Speedflow neonate IV filter 0,2 – 0,2 positive - 1,2 – 5,0 µm

Rev. 05


Baby Speedflow neonate IV filter



| | |
|-----------------------------------|---|
| PRODUCT DESCRIPTION | <p>Baby Speedflow is a non-sterile, non-toxic, self venting, 120 hours filtration (0,2 positive) / 96 hour filtration (0,2 micron) / 24 hour filtration (1.2 micron), single use device with hydrophilic PES membrane (0.2 , 0.2 positive, 1.2 or 5.0 µm) and hydrophobic PTFE membrane (0,03 µm) in a MBS housing.</p> <p>The product is provided in bulk packs for further manufacturing, processing, or repackaging.</p> |
| INTENDED USE / APPLICATION | <p>The filter is designed for use in filtration of intravenous or other aqueous solutions for removal of particles larger than 0.2 µm / 1,2 µm / 5.0 µm.</p> |
| MATERIALS | <p>Filter media: Hydrophilic PES membrane 0.2 µm / positive 0.2 µm / 1.2 µm / 5,0 µm Vent: Hydrophobic PTFE 0.03 µm Housing: Clear Acrylic-based multipolymer compound, amber or blue masterbatch.</p> <p>Inlet/Outlet connectors: Microbore tubing + double luer lock RS051 – ID 2.0mm RS052 – ID 2.2mm RS053 – ID 2.3mm RS054 – ID 2.4mm RS055 – ID 2.5mm RS056 – ID 2.8mm RS057 – ID 3.0mm RS062 – ID 2.85mm RS078 – ID 3.175mm RS079 – ID 2,7mm</p> <p>RS058 – Female Luer Lock inlet / Male Rotating Luer Lock outlet in compliance with ISO80369-7</p> <p>Flow direction: see section “Instruction / Warning”</p> |
| PRODUCT CHARACTERISTICS | <p>Dimensions WxLxD: 15.3x21.9x4.0 mm (filter body) Weight 1.35 gr. (1.7 gr. for double LL version) Hydrophilic filtration area 1.45 cm² Hydrophobic filtration area 0.25 cm² Air Flow Rate ~ 20 scc/min @ 100mbar (hydrophobic membrane) Max operating pressure 5.2 bar (75.4 psi) Max operating temperature 55 °C (131 °F)</p> <p>Minimum Water Bubble Point: PES 0.2/0.2pos µm: 3.7÷4,8 bar PES 1.2 µm: 0.7 ÷1,0 bar PES 5.0 µm: 0.15÷ 0,3 bar</p> <p>Minimum Water Flow Rate: PES 0.2pos µm : ≥3,5 ml/min @ 80 cm (31.5 in) water head pressure PES 0.2 µm : ≥4 ml/min @ 80 cm (31.5 in) water head pressure PES 1.2 µm : ≥30 ml/min @ 80 cm (31.5 in) water head pressure PES 5.0 µm : ≥55 ml/min @ 80 cm (31.5 in) water head pressure</p> <p>Bacterial Retention Brevundimonas diminuta / Candida Albicans (PES 1.2) / Not available (PES 5.0) Priming volume < 0.35 ml Pyrogenicity < 0.06 EU/ml using the LAL test method Low binding test: performed with Piperacillin Sodium, Insulin, Paclitaxel, Lidocaine HCL, Nitro-glycerin, Sodium Citrate.</p> |

Product PN RS051- RS052- RS053- RS054- RS055- RS056- RS057- RS058 - RS062 – RS078 – RS079 Mod. 984 c

Description Baby Speedflow neonate IV filter 0,2 – 0,2 positive - 1,2 – 5,0 µm Rev. 05

| | |
|--|---|
| <p>INSTRUCTIONS</p> <p>WARNINGS</p> | <p>Suggestion for easy priming procedure: keep Speedflow dry and in vertical position with the flow arrow (on the two sides of the filter) upwards. The filter will eliminate air and let the liquid flow go through. After priming is complete Speedflow filter can stay any position.</p> <p>Filter for medical use, to be assembled in clean room. Remove the external bag before planting into a clean room. Handle with care.</p> <p>Cyclohexanone for glueing is recommended. Nevertheless, if PES hydrophilic membrane comes in contact with it, membrane breaks down.</p> <p>Verify compatibility of drugs to use with the raw materials declared in specifications.</p> <p>It is not recommended to use any kind of disinfectant in direct contact with the filter. For more details, please contact GVS.</p> <p>Usage with electric/mechanical pumps - When using Speedflow Filters with any pump model, always arrange pump section above the filter and preferably keep at least 50cm between pump section and filter inlet connector.</p> <div style="text-align: center;"> <p><u>Flow direction</u></p>  </div> |
| <p>STERILIZATION</p> | <p>Ethylene oxide (Max 55°C) and gamma irradiation (Max 25 kGy).</p> |
| <p>APPLICABLE STANDARDS AND REGULATIONS</p> | <p>FOR RAW MATERIALS USED TO PRODUCE COMPONENTS:</p> <p>Test performed in compliance with USP class VI and/or ISO 10993-1.</p> <p>All materials are DEHP free, Latex free and BSE/TSE free</p> <p>Chemical composition complies with the recommendation or regulation for food contact applications.</p> <p>USA - Code of Federal Regulations, issued by Food and Drug Administration (FDA) paragraph 21 CFR 177.1500 (nylon resins).</p> <p>Test report available at GVS premises.</p> |
| <p>PACKAGING AND LABELLING</p> | <p>Box of 2.000 pcs. 2 inner PE bags of 1.000 pcs. each, Bags are separately hot sealed.</p> <p>3 bags per box (6.000 units per box).</p> <p>The first bar-code label is outside the 2 bags.</p> <p>The second bar-code label is stuck outside the box.</p> <p>Each bag is labeled with the following traceability information:</p> <ul style="list-style-type: none"> - Quantity - Product description - Product date - Lot number (OL and 5 digit batch number to trace back to raw materials used) - Operator code <p>Different lot of goods in one shipments are packed in a manner to prevent mix-ups.</p> <p>Different lot in one box are separately closed and separately labeled to prevent mix-ups.</p> |
| <p>CERTIFICATE OF COMPLIANCE</p> | <p>Conformity declaration is printed on every invoice and Certificate is according to UNI EN 10204 type 2.1</p> <p>The Quality management system is in compliance with ISO 9001:2000, ISO 13485:2003, ISO/TS 16949</p> |

Product PN RS051- RS052- RS053- RS054- RS055- RS056- RS057- RS058 - RS062 – RS078 – RS079

Mod. 984 c

Description Baby Speedflow neonate IV filter 0,2 – 0,2 positive - 1,2 – 5,0 µm

Rev. 05

| DRAWING | The attached drawing is part of this material specification and must not be duplicated or made accessible to a third party without prior written GVS SpA consent. | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
|---------------------------------|--|-----|---|-------------------------|--|-----|---------------|---|---|-----|---|---|---|-----|--------------------|---|---|-----|---|-----|---|---|-----|---|--|-----|----------------------|---|--|---|-----------------------------------|-----|---|----------------------------------|-----|----|-----------------------------|-----|----|--|-----|----|--|-----|----|--|-----|----|---|-----|----|---|-----|----|--|-----|
| VISUAL REQUIREMENTS | <p>Visual acceptance requirements apply when inspected under below conditions:</p> <p>Magnification: unaided eye, approximately 45 cm (18”) from eye Illumination: 1000 ± 200 lx or equivalent Light type: Fluorescent</p> <table border="1" data-bbox="336 640 1497 1386"> <thead> <tr> <th colspan="2">Acceptance Requirements</th> <th>AQL</th> <th>Sampling Plan</th> </tr> </thead> <tbody> <tr><td>1</td><td>Incomplete plastic support (functional)</td><td>0.1</td><td rowspan="16">ISO 2859 part. 1 1st Level</td></tr> <tr><td>2</td><td>Incomplete or misplaced membrane</td><td>0,1</td></tr> <tr><td>3</td><td>Incomplete plastic support (not functional)</td><td>0.4</td></tr> <tr><td>4</td><td>Damages, cracks or deformation on the pieces (functional)</td><td>0,1</td></tr> <tr><td>5</td><td>Damages, cracks or deformation on the pieces (non functional)</td><td>0,4</td></tr> <tr><td>6</td><td>Foreign material / Contamination > 0.2 mm²</td><td>0,1</td></tr> <tr><td>7</td><td>Embedded particles < 0.2 mm² * (max 3 per viewing area) - TAPPI DIRT ESTIMATION CHART</td><td>0,4</td></tr> <tr><td>8</td><td>Air bubbles > 0.7 mm²</td><td>0,4</td></tr> <tr><td>9</td><td>Fitting / Burr at the connection</td><td>0,4</td></tr> <tr><td>10</td><td>Burrs > 1,0 mm²</td><td>0,1</td></tr> <tr><td>11</td><td>Projecting threads from external and cones (burrs)</td><td>0.4</td></tr> <tr><td>12</td><td>Dents leaving traces, porosity, scratches.</td><td>0.4</td></tr> <tr><td>13</td><td>Plastics residual or internal membrane threads</td><td>0.4</td></tr> <tr><td>14</td><td>No loose foreign particulate upstream of the filter, plastics particles or internal membrane threads (upstream)</td><td>0.4</td></tr> <tr><td>15</td><td>Incomplete printing - pore size not readable (functional)</td><td>0,1</td></tr> <tr><td>16</td><td>Printing with smudges (max 3 < 0,2 mm² or max 5 < 0,05 mm²) - TAPPI DIRT ESTIMATION CHART</td><td>0.4</td></tr> </tbody> </table> <p>Timings: 5 sec per unit (*) Embedded Particulate Matter: according to Dirt Estimation Chart (Tappi Standard). Contamination Loose PM: free of visible particles > 0,2 mm²</p> | | | Acceptance Requirements | | AQL | Sampling Plan | 1 | Incomplete plastic support (functional) | 0.1 | ISO 2859 part. 1 1 st Level | 2 | Incomplete or misplaced membrane | 0,1 | 3 | Incomplete plastic support (not functional) | 0.4 | 4 | Damages, cracks or deformation on the pieces (functional) | 0,1 | 5 | Damages, cracks or deformation on the pieces (non functional) | 0,4 | 6 | Foreign material / Contamination > 0.2 mm ² | 0,1 | 7 | Embedded particles < 0.2 mm ² * (max 3 per viewing area) - TAPPI DIRT ESTIMATION CHART | 0,4 | 8 | Air bubbles > 0.7 mm ² | 0,4 | 9 | Fitting / Burr at the connection | 0,4 | 10 | Burrs > 1,0 mm ² | 0,1 | 11 | Projecting threads from external and cones (burrs) | 0.4 | 12 | Dents leaving traces, porosity, scratches. | 0.4 | 13 | Plastics residual or internal membrane threads | 0.4 | 14 | No loose foreign particulate upstream of the filter, plastics particles or internal membrane threads (upstream) | 0.4 | 15 | Incomplete printing - pore size not readable (functional) | 0,1 | 16 | Printing with smudges (max 3 < 0,2 mm ² or max 5 < 0,05 mm ²) - TAPPI DIRT ESTIMATION CHART | 0.4 |
| Acceptance Requirements | | AQL | Sampling Plan | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 1 | Incomplete plastic support (functional) | 0.1 | ISO 2859 part. 1 1 st Level | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 2 | Incomplete or misplaced membrane | 0,1 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 3 | Incomplete plastic support (not functional) | 0.4 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 4 | Damages, cracks or deformation on the pieces (functional) | 0,1 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 5 | Damages, cracks or deformation on the pieces (non functional) | 0,4 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 6 | Foreign material / Contamination > 0.2 mm ² | 0,1 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 7 | Embedded particles < 0.2 mm ² * (max 3 per viewing area) - TAPPI DIRT ESTIMATION CHART | 0,4 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 8 | Air bubbles > 0.7 mm ² | 0,4 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 9 | Fitting / Burr at the connection | 0,4 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 10 | Burrs > 1,0 mm ² | 0,1 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 11 | Projecting threads from external and cones (burrs) | 0.4 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 12 | Dents leaving traces, porosity, scratches. | 0.4 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 13 | Plastics residual or internal membrane threads | 0.4 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 14 | No loose foreign particulate upstream of the filter, plastics particles or internal membrane threads (upstream) | 0.4 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 15 | Incomplete printing - pore size not readable (functional) | 0,1 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 16 | Printing with smudges (max 3 < 0,2 mm ² or max 5 < 0,05 mm ²) - TAPPI DIRT ESTIMATION CHART | 0.4 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| PERFORMANCE REQUIREMENTS | <table border="1" data-bbox="349 1509 1477 1948"> <thead> <tr> <th colspan="2">Acceptance Requirement</th> <th>AQL</th> <th>Sampling Plan</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>Bubble point to verify PES integrity</td> <td rowspan="3">0,1</td> <td rowspan="3">ISO 2859 part. 1 1st Level</td> </tr> <tr> <td></td> <td>- 0.2 / 0.2pos µm: 3.7 ÷ 4,8 bar (ramped pressure in 15 seconds) - 1.2 µm: 0.7 ÷ 1,0 bar - 5.0 µm: 0.15 ÷ 0,3 bar</td> </tr> <tr> <td>2</td> <td>WBT to verify PTFE</td> </tr> <tr> <td>3</td> <td>Burst test to verify housing pressure integrity</td> <td>0,1</td> <td></td> </tr> <tr> <td>4</td> <td>Water Flow rate @ 80 cm water head pressure</td> <td rowspan="4">0,1</td> <td rowspan="4"></td> </tr> <tr> <td></td> <td>- 0.2pos µm: ≥ 3,5 ml/min</td> </tr> <tr> <td></td> <td>- 0.2 µm: ≥ 4 ml/min</td> </tr> <tr> <td></td> <td>- 1.2 µm: ≥ 30 ml/min - 5.0 µm: ≥ 55 ml/min</td> </tr> </tbody> </table> | | | Acceptance Requirement | | AQL | Sampling Plan | 1 | Bubble point to verify PES integrity | 0,1 | ISO 2859 part. 1 1 st Level | | - 0.2 / 0.2pos µm: 3.7 ÷ 4,8 bar (ramped pressure in 15 seconds) - 1.2 µm: 0.7 ÷ 1,0 bar - 5.0 µm: 0.15 ÷ 0,3 bar | 2 | WBT to verify PTFE | 3 | Burst test to verify housing pressure integrity | 0,1 | | 4 | Water Flow rate @ 80 cm water head pressure | 0,1 | | | - 0.2pos µm: ≥ 3,5 ml/min | | - 0.2 µm: ≥ 4 ml/min | | - 1.2 µm: ≥ 30 ml/min - 5.0 µm: ≥ 55 ml/min | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Acceptance Requirement | | AQL | Sampling Plan | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 1 | Bubble point to verify PES integrity | 0,1 | ISO 2859 part. 1 1 st Level | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | - 0.2 / 0.2pos µm: 3.7 ÷ 4,8 bar (ramped pressure in 15 seconds) - 1.2 µm: 0.7 ÷ 1,0 bar - 5.0 µm: 0.15 ÷ 0,3 bar | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 2 | WBT to verify PTFE | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 3 | Burst test to verify housing pressure integrity | 0,1 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| 4 | Water Flow rate @ 80 cm water head pressure | 0,1 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | - 0.2pos µm: ≥ 3,5 ml/min | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | - 0.2 µm: ≥ 4 ml/min | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | - 1.2 µm: ≥ 30 ml/min - 5.0 µm: ≥ 55 ml/min | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |

This material specification describes the properties of product above indicated.
This document contains general requirements, material description, drawing references, defect specification, biological material requirements.

PRODUCT SPECIFICATION

| | | |
|-------------|--|------------|
| Product PN | RS051- RS052- RS053- RS054- RS055- RS056- RS057- RS058 - RS062 – RS078 – RS079 | Mod. 984 c |
| Description | Baby Speedflow neonate IV filter 0,2 – 0,2 positive - 1,2 – 5,0 µm | Rev. 05 |

REVISIONS AND APPROVALS:

| DATE | REV. | REASON FOR CHANGE | ISSUED AND CONTROLLED BY: (name /function and signature) | APPROVED BY: (name /function and signature) |
|------------|------|--|---|---|
| 14/06/2022 | 12 | Updated limit for LAL test from < 0.25 EU/ml to < 0.06 EU/ml. Added code. RS056BCYRH050M00 for Baby Speedflow amber 5 micron. The amber version is applicable to all versions / pore size. | Elsa Caruso - RP <i>Elsa Caruso</i> | Barbara Finessi - AQP <i>Barbara Finessi</i> Enrico Salvarani – RPROG <i>Enrico Salvarani</i> Luca Zanini – DAM HC <i>Luca Zanini</i> Tiziana Landi – DAQ <i>Tiziana Landi</i> |

Customer Approval:

We accept this material specification as a part of the agreed terms of delivery

Company name _____

Approved by: _____
(Name, Function) (Signature)

Date _____
(Company stamp)

Please send back this document signed for approval. If we will not receive this specification signed, we consider the first order placed as implicit approval.