ENGLISH Instructions for Use

2800/R1, 2800/R2 and 2800/R3 Re-usable Filters

Intended use: Spirometer pre filters are used in pulmonary function tests for trapping bacteria and viruses and other particulates, reducing the risk of cross-contamination between the patient and the machine.

Bidirectional Filter.

Code for replacement pads: 7270/007.

Place clean pad in the moulding case, so the pad sits inside the internal housing.

Place the lid moulding on top of the moulding case, so the internal seal is located.

Apply pressure and pinch both mouldings together.

Cautions

CAUTION: A CAUTION statement provides important information about a potentially hazardous situation which, if not avoided, may result in harm to the patient.

Any use of the medical device requires full understanding and strict observation of all portions of these Instructions for Use. The medical device may only be used for the purpose specified under "Intended Use". Observe all CAUTION statements throughout this manual and all statements on medical device labels. The manufacturer disclaims any liability for patient injury caused by noncompliance with these statements.

Following installation of the filter, check that all connections are secure.

The medical device IFU is not available for each individual filter. Only one copy of the Instructions for Use is included in the distribution package and should therefore be kept in an accessible location for users.

Medical device is only to be used by healthcare professionals.

User has the responsibility to report any serious incident that has occurred in relation to the medical device, to the manufacturer, EC REP and the competent authority of the Member State in which the user and/or the patient is established.

To avoid contamination and soiling, the product should remain packaged until ready to be used.

Do not use the product if the packaging and/or product is damaged.

If product was kept outside the stated storage temperature, product functionality cannot be guaranteed. Storage Temperature: +5°C to 40°C.

Store the product at room temperature, in dry place, free of dust and not exposed to direct sunlight.

Handle filter with care. Abusive handling can damage the filter.

For Your Safety and that of Your Patients, follow the Instruction for Use of the machine!

Before installation, check that all system components are free of obstructions and of foreign bodies. Otherwise, use is limited, or faulty operation is possible.

Inappropriate connection to other devices can result in a patient hazard.

Do not use after expiry date stated on primary packaging.

Manufacturing date

Expiry date

Monitor filter performance during the use. If any changes to filter performance or appearance are observed, change the pad to new one.

Always keep spare pad ready to exchange with the one in use, if any of the below are observed:

- Increase in resistance;
- Filter is visibly soiled or damaged or if it rattles;
- Excessive condensation;
- Associated alarm.

Single use of the pad.

(Symbol not applicable for re-usable plastic parts).

The pad must not be reused. Reuse, cleaning or sterilization may lead to failure of the products, risk infection and injury to the patient.

The device may only be cleaned by the user in accordance to the below manufacturer's instructions.

Cleaning instructions:

1. Take the assembled product apart and discard the filtration pad in accordance with the local hospital hygiene and waste disposal regulations for hazardous infectious clinical waste.
2. Hand wash the plastic parts in hot water with mild soap.
3. Rinse all parts in clean water, shake off excess water and allow to air dry.
4. Wipe parts with alcohol impregnated wipes.
5. Plastic parts must be cleaned after each patient.

This product is latex and PHT free.

Batch code symbol.

Prescription only.

Manufacturer.

Importer.

Distributor.

Operator's manual. Attention: Consult instructions for use.

Unique Device Identifier.

Packaging unit. To indicate the number of pieces inside the package.

Authorized representative in the European Union.

Following use, the products must be disposed off in accordance with the local hospital hygiene and waste disposal regulations for hazardous infectious clinical waste.

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