GVS Filter Technology UK Ltd., Caton Road, Lancaster, LA1 3PE, United Kingdom, tel. +44 (0)1524 847600 www.gvsuk@gvs.com

English Instructions for Use

2802 series

Intended use: Spiroguard Machine Adaptor acts as an adaptor between the spiroguard filter and the pulmonary function machine.

No limitation to the product user population.

The device may only be cleaned by the user in accordance to the below manufacturer's instructions.

Hand wash adaptor with hot water and a mild soap. Shake it gently to remove trapped water droplets. Air dry completely and wipe with alcohol impregnated wipes before reconnecting it to the machine. **Adaptor must be cleaned after each patient.**

Cautions

CAUTION: A CAUTION statement provides important information about a potentially hazardous situation which, if not avoided, may result in harm to the patient.

Any use of the medical device requires full understanding and strict observation of all portions of these Instructions for Use. The medical device may only be used for the purpose specified under "Intended Use". Observe all CAUTION statements throughout this manual and all statements on medical device labels. The manufacturer disclaims any liability for patient injury caused by noncompliance with these statements.

Following installation of the filter, check that all connections are secure.

The medical device IFU is not available for each individual adaptor. Only one copy of the Instructions for Use is included in the distribution package and should therefore be kept in an accessible location for users.


Medical device is only to be used by healthcare professionals.

User has the responsibility to report any serious incident that has occurred in relation to the medical device, to the manufacturer, EC REP and the competent authority of the Member State in which the user and/or the patient is established.

To avoid contamination and soiling, the product should remain packaged until ready to be used.

Do not use the product if the packaging and/or product is damaged.

If product was kept outside the stated storage temperature, product functionality cannot be guaranteed. Storage Temperature: +5°C to 40°C.

Store the product at room temperature, in dry place, free of dust and not exposed to direct sunlight.

Handle product with care. Abusive handling can damage the filter.

For Your Safety and that of Your Patients, follow the Instruction for Use of the machine!

Before installation, check that all system components are free of obstructions and of foreign bodies. Otherwise, use is limited, or faulty operation is possible.

Inappropriate connection to other devices can result in a patient hazard.

Do not use after expiry date stated on primary packaging.

Manufacturing date

Expiry date

This product is latex and PHT free.

Batch code symbol.

Prescription only.

Manufacturer.

Importer.

Distributor.

Operator's manual. Attention: Consult instructions for use.

Unique Device Identifier.

Packaging unit. To indicate the number of pieces inside the package.

Authorized representative in the European Union.

Following use, the products must be disposed off in accordance with the local hospital hygiene and waste disposal regulations for hazardous infectious clinical waste.

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