




FILTER TECHNOLOGY

BACTERIOLOGICAL PETRI DISHES

Petri Dishes, Plastic Material

- ◆ Non-cytotoxic virgin polystyrene, excellent optical clarity
- ◆ Engineered for optimum flatness to provide uniform agar thickness
- ◆ Single or more vents are available: no-vent for anaerobic and long-term work, single vent for limiting gas exchange, multi vents for improving gas exchange



General petri Dishes

90 and 150mm Petri dishes are manufactured in accordance with standard ISO 24998:2008, ideal for use with most automatic filling machine.

Product Code	Description	Quantity
PETRIDPC035151NVS	Ø 35x15 mm, 1 Room	10 pcs / polybag, 300 pcs / carton
PETRIDPC049131NVS	Ø 49 x 13 mm, 1 Room	10 pcs / polybag, 1,000 pcs / carton
PETRIDPC055151NVS	Ø 55x15 mm, 1 Room	10 pcs / polybag, 1,000 pcs / carton
PETRIDPC0601513VS	Ø 60x15mm, 1 Room / 3 vents	10 pcs / polybag, 1,000 pcs / carton
PETRIDPC070151NVS	Ø 70 x 15 mm, 1 Room	20 pcs / polybag, 800 pcs / carton
PETRIDPC0901513VS	Ø 90 x 15 mm, 1 Room / 3 vents	10 pcs / polybag, 500 pcs / carton
PETRIDPC1501513VS	Ø 150 x 15 mm, 1 Room / 3 vents	10 pcs / polybag, 200 pcs / carton

Compartmented Dishes

Ideal for different media on same dish, or meet the requirement of saving incubator space.

Product Code	Description	Quantity
PETRIDPC0901523VS	Ø 90 x 15mm, 2 Rooms / 3 vents	10 pcs / polybag, 500 pcs / carton
PETRIDPC0901533VS	Ø 90 x 15mm, 3 Rooms / 3 vents	10 pcs / polybag, 500 pcs / carton

Square Dishes

Noncompartmented, ideal for antibiotic sensitivity testing.

Product Code	Description	Quantity
PETRIDPS130151NVS	Ø 130 x 130 x 15 mm, 1 Room	10 pcs / polybag, 150 pcs / carton

Contact Plates

10x10mm counting net, alphanumeric, facilitates colony counting and locating.

Product Code	Description	Quantity
PETRIDPC0651513VS	Ø 65 x 15 mm, 1 Room / 3 vents	10 pcs / polybag, 1,000 pcs / carton

Petri Dishes, Glass Material

3.3 Borosilicate glass, withstand repeated sterilization (wet or dry).



Product Code	Description	Quantity
PETRIDGC060151NVS	Ø 60 mm x 15 mm	288/pk
PETRIDGC075151NVS	Ø 75 mm x 15 mm	288/pk
PETRIDGC090181NVS	Ø 90 mm x 18 mm	192/pk
PETRIDGC100201NVS	Ø 100 mm x 20 mm	192/pk

Petri Dish Carrier Rack

Designed for holding and carrying plastic or glass petri dishes.

- ◆ Made of STAINLESS STEEL Material
- ◆ Autoclavable and resistant to universal chemicals
- ◆ Highly polished surface, no edges and sharps
- ◆ Lightweight, easy to load, safe for transporting

Packing: 1 pc / inner box, 5 pcs / carton

Product Code	Dimensions	Description
PETRIDC340110380330	340 × 110 × 380 mm	3 sections, holds 30 petri dishes

