

Product PN 10486001, 10486081

Mod. 984

Description FAST ACCESSORIES

Rev. 02

FAST ACCESSORIES – Slide Holders



PRODUCT DESCRIPTION

FAST Slide Accessories for protein array processing.

| FAST Slide Holders | Item Number |
|-------------------------|-------------|
| FAST Frame Slide Holder | 10486001 |
| Chip Clip Slide Holder | 10486081 |

MANUFACTURER NAME

GVS North America, Inc.
63 Community Drive
Sanford, ME 04073
Phone: +1.866.736.1250
Email: CustomerCareGVSLs@gvs.com
Website: www.gvs.com

INTENDED USE / APPLICATION

- FAST Slide Holders
 - FAST Frame holds up to 4 FAST slides (with chambers) for protein microarray processing.
 - Chip Clip hold one FAST slide (with chamber) for processing microarrays.

PRODUCT CHARACTERISTIC

Regulatory Documentation Supplied:

- ☐ Biocompatibility according ISO 10993-1
- ☐ IMDS
- ☒ RoHS, directive 2002/95/CE
- ☒ BSE/TSE, directive 2003/32/CE
- ☐ DEHP plasticizer Free and latex free
- ☐ Aging
- ☒ Reach 1907/2006/CE (hazardous substances regulation)
- ☐ Dir. 67/548/CE and Reg. 1272/2008/CE (medical sector dangerous substances)
- ☒ Conflict minerals

| Fast Slide Holder Number | Number of Slides | Footprint (mm) |
|--------------------------|------------------|----------------|
| 10486001 | Up to 4 | 128 x 86 |
| 10486081 | one | 50 x 85 |


PRODUCT SPECIFICATION

Product PN 10486001, 10486081

Mod. 984

Description FAST ACCESSORIES



Rev. 02

| | |
|--------------------------------|---|
| MATERIALS | <ul style="list-style-type: none"> FAST Slide Holders <ul style="list-style-type: none"> Made of autoclavable plastic, compatible with standard 1" x 3" glass slides |
| STORAGE | Store at recommended room temperature. |
| PRODUCT SHELF LIFE | Slide holders: Use is unlimited |
| PACKAGING AND LABELLING | <p>Packaging: Primary: FAST slide holders are placed in a labeled box. Example of label</p>  |

This material specification describes the properties of product above indicated.

This document contains general requirements, material description, drawing references, defect specification, biological material requirements.

REVISIONS AND APPROVALS:

| DATE | REV. | REASON FOR CHANGE | ISSUED AND CONTROLLED BY: (name /function and signature) | APPROVED BY: (name /function and signature) |
|--------|------|-------------------|--|--|
| 6/1/18 | 0 | Initial release | Debra English/Manager Biochemistry Applications  | Kevin Wrigley/Director, Quality  |

Customer Approval:

PRODUCT SPECIFICATION

Product PN 10486001, 10486081

Mod. 984

Description FAST ACCESSORIES

Rev. 02

We accept this material specification as a part of the agreed terms of delivery

Company name _____

Approved by: _____
(Name, Function) (Signature)

Date _____
(Company stamp)

Please send back this document signed for approval. If we will not receive this specification signed, we consider the first order placed as implicit approval.