

Product PN 10485331


Mod. 984

Description

Protein Arraying Buffer (2X); 4/pk

Rev. 02

Protein Arraying Buffer (2X); 4/pk



PRODUCT DESCRIPTION	Protein array buffers are optimized for use on FAST Slides. Protein Arraying Buffer is supplied as a 2X concentrate in 10 ml plastic bottles.
MANUFACTURER NAME	GVS North America, Inc. 63 Community Drive Sanford, ME 04073 Phone: +1.866.736.1250 eMail: CustomerCareGVSLs@gvs.com - Website: www.gvs.com
INTENDED USE / APPLICATION	Protein Arraying Buffer enhances long-term protein stability and molecular recognition activity of arrayed proteins. Protein Arraying Buffer is a 2X concentrate which is diluted to 1X arraying solution by the addition of protein.
MATERIALS	<p>Liquid: Protein Arraying Buffer</p> <p>Color: 2 X Concentrate Clear liquid</p> <p>Bottle: polypropylene</p> <p>Regulatory Documentation Required:</p> <ul style="list-style-type: none"> <input type="checkbox"/> Biocompatibility according ISO 10993-1 <input type="checkbox"/> IMDS <input checked="" type="checkbox"/> Rohs, directive 2002/95/CE <input checked="" type="checkbox"/> BSE/TSE, directive 2003/32/CE <input type="checkbox"/> DEHP plasticizer Free and latex free <input type="checkbox"/> Aging <input checked="" type="checkbox"/> Reach 1907/2006/CE (hazardous substances regulation) <input type="checkbox"/> Dir. 67/548/CE and Reg. 1272/2008/CE (medical sector dangerous substances) <input type="checkbox"/> Conflict minerals
PRODUCT CHARACTERISTIC	<p>2 X Concentrate Clear liquid dispensed in polypropylene bottle.</p> <p>Volume/bottle: 10 ml pH: 7.4 ± 0.1</p> <p>Toxicity: See SDS for further information.</p> <p>Storage Temperature: Room temperature. Ensure cap is on tight.</p>

PRODUCT SPECIFICATION


Product PN 10485331

Mod. 984

Description

Protein Arraying Buffer (2X); 4/pk

Rev. 02

<p>INSTRUCTIONS WARNINGS</p>	<p>Contains methanol. See SDS for further information.</p>														
<p>PRODUCT SHELF LIFE</p>	<p>One year</p>														
<p>PACKAGING AND LABELLING</p>	<p>Packaging: Primary: Protein Arraying Buffer dispensed in PET bottle. Label placed on bottle. Secondary: Bottle placed in labeled box.</p> 														
<p>CERTIFICATE OF COMPLIANCE</p>	<p>The Quality management system is in compliance with ISO 9001. When applicable, GVS NA will provide the CoC with item number, lot number and expiration date.</p>														
<p>VISUAL REQUIREMENTS</p>	<p>Visual acceptance requirements apply when inspected under below conditions: Magnification: none Light type: <i>normal room lighting</i></p> <table border="1" data-bbox="363 1272 1458 1355"> <thead> <tr> <th colspan="2">Acceptance Requirement</th> <th>AQL</th> <th>Sampling Plan</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>Clear liquid – no precipitate</td> <td>All</td> <td>100%</td> </tr> </tbody> </table>			Acceptance Requirement		AQL	Sampling Plan	1	Clear liquid – no precipitate	All	100%				
Acceptance Requirement		AQL	Sampling Plan												
1	Clear liquid – no precipitate	All	100%												
<p>PERFORMANCE REQUIREMENTS</p>	<table border="1" data-bbox="437 1442 1230 1606"> <thead> <tr> <th colspan="3">Acceptance Requirement</th> </tr> <tr> <th></th> <th>Parameter</th> <th>Specification</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>pH</td> <td>7.4 ± 0.1</td> </tr> <tr> <td>2</td> <td>Conductivity</td> <td>10 – 18.4 mS</td> </tr> </tbody> </table>			Acceptance Requirement				Parameter	Specification	1	pH	7.4 ± 0.1	2	Conductivity	10 – 18.4 mS
Acceptance Requirement															
	Parameter	Specification													
1	pH	7.4 ± 0.1													
2	Conductivity	10 – 18.4 mS													

This material specification describes the properties of product above indicated.
This document contains general requirements, material description, drawing references, defect specification, biological material requirements.



PRODUCT SPECIFICATION

Product PN 10485331


Mod. 984

Description

Protein Arraying Buffer (2X); 4/pk

Rev. 02

REVISIONS AND APPROVALS:

DATE	REV.	REASON FOR CHANGE	ISSUED AND CONTROLLED BY: (name /function and signature)	APPROVED BY: (name /function and signature)
11/20/2019	0	Initial Release	Debra English/Manager Biochemistry Applications 	Kevin M. Wrigley, Director, Quality 

Customer Approval:

We accept this material specification as a part of the agreed terms of delivery

Company name _____

Approved by: _____
(Name, Function) (Signature)

Date _____
(Company stamp)

Please send back this document signed for approval. If we will not receive this specification signed, we consider the first order placed as implicit approval.