


Product PN	1227204, 1227205, 1229451, 1226150, 1227207, 1227208, 1227209	Mod. 984
Description	30 mm Cameo Syringe Filters w/Glass Fiber Membrane	Rev. 02



30 mm Cameo, Glass Fiber (GF)

PRODUCT DESCRIPTION	Non-Sterile 30mm Syringe Filter made of Polypropylene Housing assembled with various pore sizes of Glass Fiber/Polypropylene (GF/PP) Membrane							
	Membrane Material	Pore Size (um)	Color	Packaging 50/pk	Packaging 200/pk	Packaging 500/pk	Packaging 1000/pk	Packaging 5000/pk
	Glass Fiber (GF)	1.00	Transparent	1227204	N/A	1227205	1229451	1226150
	Glass Fiber (GF)	0.70	Transparent	1227207	N/A	N/A	1227208	1227209
MANUFACTURER NAME	GVS North America 63 Community Drive Sanford, Me 04073 Phone: +1.866.736.1250 eMail: CustomerCareGVSLs@gvs.com - Website: www.gvs.com							
INTENDED USE / APPLICATION	Applications <ul style="list-style-type: none"> Analytical sample preparation Dissolution testing Content uniformity Environmental samples Composite assays Food analysis Biofuel analysis 							
MATERIALS	Filter media: Glass Fiber Frame/Housing Polymer: Polypropylene Color: Transparent Other insert(s): N/A Regulatory Documentation Required: <ul style="list-style-type: none"> <input type="checkbox"/> Biocompatibility according ISO 10993-1 <input type="checkbox"/> IMDS <input type="checkbox"/> DEHP plasticizer Free and latex free <input checked="" type="checkbox"/> Rohs, Directive 2002/32/CE <input type="checkbox"/> Aging <input checked="" type="checkbox"/> BSE/TSE, directive 2003/32/CE <input checked="" type="checkbox"/> 1907/2006/CE (hazardous substances regulation) <input type="checkbox"/> Dir. 67/548/CE and Reg. 1272/2008/CE (medical sector dangerous substances) <input checked="" type="checkbox"/> Conflict minerals 							

Product PN	1227204, 1227205, 1229451, 1226150, 1227207, 1227208, 1227209	Mod. 984
Description	30 mm Cameo Syringe Filters w/Glass Fiber Membrane	Rev. 02

PRODUCT CHARACTERISTIC	<p>Characteristics</p> <p>Membrane Material: Glass Fiber Membrane Diameter: 30 mm Effective Filtration Area: 4.8 cm² Housing Diameter: 33 mm Housing Material: Pure polypropylene is heat-sealed without the use of glues or sealants Inlet / Outlet: FLL-MLS Prefilter: 1.0 μ binderless glass-fiber, in some configurations Holdup Volume: <120 microliter Maximum Operating Temperature: 82°C/180°F Maximum Operating Pressure: 80 psi Sterile: No</p>
PRODUCT SHELF LIFE	When stored under normal storage conditions, this product should be stable for 5 years
STERILIZATION	<p>This is a non-sterile product, compatible with these sterilization techniques*:</p> <ul style="list-style-type: none"> <input checked="" type="checkbox"/> EtO <input checked="" type="checkbox"/> Gamma (up to 40 kGy) <input type="checkbox"/> Autoclave (121°C) <input type="checkbox"/> eBeam <p><i>*user has the responsibility to validate their process</i></p>
COMPLIANCE	The Quality management system is in compliance with ISO 9001:2015
DRAWING	 <p>1. FEMALE LUER LOCK INLET 2. MEMBRANE 3. MALE LUER SLIP OUTLET 4. OVERMOLDED RING</p> <p>LUER CONNECTIONS CONFORM TO ISO 594</p> <p>Ø33.045 22.860</p> <p>DIMENSIONS ARE IN MM, FOR REFERENCE ONLY</p>



PRODUCT SPECIFICATION

Product PN	1227204, 1227205, 1229451, 1226150, 1227207, 1227208, 1227209	Mod. 984
Description	30 mm Cameo Syringe Filters w/Glass Fiber Membrane	Rev. 02

VISUAL REQUIREMENTS	<i>Visual acceptance requirements apply when inspected under below conditions:</i>																																																										
	Magnification:																																																										
	Light type:																																																										
	Timings:																																																										
	<table border="1"> <thead> <tr> <th colspan="3">Acceptance Requirement</th> <th>Sampling Plan</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>Mold Flash</td> <td><0.2 mm</td> <td>100%</td> </tr> <tr> <td>2</td> <td>Voids/Sinks</td> <td>None</td> <td>100%</td> </tr> <tr> <td>3</td> <td>Short shots</td> <td>None</td> <td>100%</td> </tr> <tr> <td>4</td> <td>Missing/Jutting Media/Mesh</td> <td>None</td> <td>100%</td> </tr> <tr> <td>5</td> <td>Scuffing</td> <td>Total length of scuff exceed 2 ribs</td> <td>100%</td> </tr> <tr> <td>6</td> <td>Luer Damage</td> <td>None</td> <td>100%</td> </tr> <tr> <td>7</td> <td>Weld Damage</td> <td>None</td> <td>100%</td> </tr> <tr> <td>8</td> <td>Heat Stake Issues</td> <td>None</td> <td>100%</td> </tr> <tr> <td>9</td> <td>Spots Molded in Plastic</td> <td>None</td> <td>100%</td> </tr> <tr> <td>10</td> <td>Discoloring</td> <td>None</td> <td>100%</td> </tr> <tr> <td>11</td> <td>Mesh out of Place/Media Off-Center</td> <td>None</td> <td>100%</td> </tr> <tr> <td>12</td> <td>Pad Printing Incomplete</td> <td>None</td> <td>100%</td> </tr> <tr> <td>13</td> <td>Embedded Particles</td> <td>< 0.8 mm² (Maximum 3 particles)</td> <td>100%</td> </tr> </tbody> </table>			Acceptance Requirement			Sampling Plan	1	Mold Flash	<0.2 mm	100%	2	Voids/Sinks	None	100%	3	Short shots	None	100%	4	Missing/Jutting Media/Mesh	None	100%	5	Scuffing	Total length of scuff exceed 2 ribs	100%	6	Luer Damage	None	100%	7	Weld Damage	None	100%	8	Heat Stake Issues	None	100%	9	Spots Molded in Plastic	None	100%	10	Discoloring	None	100%	11	Mesh out of Place/Media Off-Center	None	100%	12	Pad Printing Incomplete	None	100%	13	Embedded Particles	< 0.8 mm ² (Maximum 3 particles)	100%
	Acceptance Requirement			Sampling Plan																																																							
	1	Mold Flash	<0.2 mm	100%																																																							
	2	Voids/Sinks	None	100%																																																							
	3	Short shots	None	100%																																																							
	4	Missing/Jutting Media/Mesh	None	100%																																																							
	5	Scuffing	Total length of scuff exceed 2 ribs	100%																																																							
	6	Luer Damage	None	100%																																																							
	7	Weld Damage	None	100%																																																							
8	Heat Stake Issues	None	100%																																																								
9	Spots Molded in Plastic	None	100%																																																								
10	Discoloring	None	100%																																																								
11	Mesh out of Place/Media Off-Center	None	100%																																																								
12	Pad Printing Incomplete	None	100%																																																								
13	Embedded Particles	< 0.8 mm ² (Maximum 3 particles)	100%																																																								
PERFORMANCE REQUIREMENTS	<table border="1"> <thead> <tr> <th colspan="3">Acceptance Requirement</th> <th>Sampling Plan</th> </tr> </thead> <tbody> <tr> <td>Pore size</td> <td>1.0 um</td> <td>0.70 um</td> <td>AQL 0.1 Special inspection level S3</td> </tr> <tr> <td>Pressure</td> <td>≥ 80 PSI 10 Seconds</td> <td>≥ 80 PSI 10 Seconds</td> <td>ANSI/ASQ Standard Z1.4 - 2008</td> </tr> <tr> <td>Bubble point</td> <td>N/A</td> <td>N/A</td> <td>N/A</td> </tr> </tbody> </table>			Acceptance Requirement			Sampling Plan	Pore size	1.0 um	0.70 um	AQL 0.1 Special inspection level S3	Pressure	≥ 80 PSI 10 Seconds	≥ 80 PSI 10 Seconds	ANSI/ASQ Standard Z1.4 - 2008	Bubble point	N/A	N/A	N/A																																								
	Acceptance Requirement			Sampling Plan																																																							
	Pore size	1.0 um	0.70 um	AQL 0.1 Special inspection level S3																																																							
	Pressure	≥ 80 PSI 10 Seconds	≥ 80 PSI 10 Seconds	ANSI/ASQ Standard Z1.4 - 2008																																																							
Bubble point	N/A	N/A	N/A																																																								



This material specification describes the properties of product above indicated.
 This document contains general requirements, material description, drawing references, defect specification, biological material requirements.



PRODUCT SPECIFICATION

Product PN	1227204, 1227205, 1229451, 1226150, 1227207, 1227208, 1227209	Mod. 984
Description	30 mm Cameo Syringe Filters w/Glass Fiber Membrane	Rev. 02

REVISIONS AND APPROVALS:

DATE	REV.	REASON FOR CHANGE	ISSUED AND CONTROLLED BY: (name /function and signature)	APPROVED BY: (name /function and signature)
11/18/19	0	Initial Release	Joe DeSisto, Director, Process Engineering 	Kevin Wrigley, Director, Quality 

Customer Approval:

We accept this material specification as a part of the agreed terms of delivery

Company name _____

Approved by: _____
 (Name, Function) (Signature)

Date _____
 (Company stamp)

Please send back this document signed for approval. If we will not receive this specification signed, we consider the first order placed as implicit approval.