


4100/735

Oxygen Concentrator HEPA Filter



Main Air Intake Refill For Sunrise 493-0012-007

Filtration Efficiency	BFE 99.999%
Casing	Clear Bronze Tinted ABS
Noise	Unique Air Intake Chamber Reduces Noise
Connections	22mm Push fit
Usage	Single Use Only

*Independently tested, data available upon request



DIMENSIONS ARE APPROXIMATE



TECHNICAL DATA SHEET

MAF-TDS-154

Issue 4
Issue Date: 27.8.13

4100/735