

PRODUCT SPECIFICATION

Glass Fiber Filters w/Binder, 1.0um



PRODUCT DESCRIPTION	Glass Fiber membranes with binders, composed of Borosilicate glass fibers			
	Part Number	Diameter	Quantity/pk	
	1215557	13 mm	100/Pk	
	1215559	25 mm	100/Pk	
	1215561	42 mm	100/Pk	
	1215562	47 mm	100/Pk	
	1215563	75 mm	25/Pk	
	1215564	90 mm	25/Pk	
	1215565	124 mm	25/Pk	
	1215566	127 mm	25/Pk	
	1215567	142 mm	25/Pk	
	1215568	257 mm	25/Pk	
	1215569	293 mm	25/Pk	
	1268603	24 x 24 cm	10/Pk	
	Sanford, ME 04073 Phone: +1.866.736.1250 Email: CustomerCareGVSLS@gvs.com Website: www.gvs.com			
INTENDED USE/ APPLICATION	Typical Applications:			
MATERIALS	Material: Quartz filter paper			
	Regulatory Documentation Re Biocompatibility according IMDS Rohs, directive 2002/95, BSE/TSE, directive 2003 DEHP plasticizer and la	ng ISO 10993-1 /CE 3/32/CE		



PRODUCT SPECIFICATION

12155		557, 1215559, 1215561, 1215562, 1215563, 1215564, 1215565, 1215566, 1215567, 568, 1215569, 1268603 5 Fiber Filters w/Binder	
PRODUCT CHARACTERIS	TICS	 Weight: 1.0-2.0 g/m² Thickness: 12 mils Acrylic binder Biologically inert Bonding reduces media migration 	
PRODUCT SHELF LIFE 5 years in normal storage conditions: Store the product at room temperate, in a dry place, free of dust and not exposed to determine the product at room temperate.		ct sunlight.	
STERILIZATION	N	Not applicable - This is a non-sterile product.	
PACKAGING AND LABELING Bulk packed filters in a gloss finished, chipboard box Each box will be labeled with: Quantity Product item code and description Lot number Storage conditions			
CERTIFICATE C	OF	A conformity declaration will be sent to the customer for each shipment with lot number a manufacture. The Quality management system is in compliance with ISO 9001.	and date of
This material specifi This document cont	ication des	scribes the properties of product above indicated. ral requirements, material description, drawing references, defect specification, and biological material	requirements.

REVISIONS AND APPROVALS:

DATE	REV	REASON FOR CHANGE	ISSUED AND CONTROLLED BY: (name /function and signature)	APPROVED BY: (name /function and signature)
11/4/19	0	Initial release	Joe DeSisto, Director, Process Engineering	Kevin Wrigley/Director of Quality
			fre Debito	Sem m. u)



PRODUCT SPECIFICATION

Product PN	1215557, 1215559, 1215561, 1215562, 1215563, 1215564, 1215565, 1215566, 1215567, 1215568, 1215569, 1268603	Mod. 984 Rev. 01
Description	Glass Fiber Filters w/Binder	Nev. 01

CUSTOMER APPROVAL:

We accept this m	aterial specification as a part of the agreed	terms of delivery.	
Company name _			_
Approved by:	(Name, Function)	(Signature)	
Date(Company stamp)			

Please send back this document signed for approval. If we will not receive this specification signed, we consider the first order placed as implicit approval.