

## PRODUCT SPECIFICATION

Product PN	3032875, 3032876, 3029505, 3018505, 1239465, 1235748, 3057106	Mod. 984
Description	0.03 μm Polyethersulfone Filtration Membrane Discs and Sheets	Rev. 01

# 0.03 µm Polyethersulfone (PES) Filtration Membrane in Discs and Sheets



PRODUCT DESCRIPTION	0.03 μm PES Discs and Sheets						
	Part Number	Geometry	Dimensions	Packaging			
	3032875	Disc	13 mm	100/Pk			
	3032876	Disc	25 mm	100/Pk			
	3029505	Disc	47mm	100/Pk			
	3018505	Disc	90 mm	25/pk			
	1239465	Sheet	150 x 50 mm	5/pk			
	1235748	Sheet	200 x 200 mm	5/pk			
	3057106	Sheet	30 cm x3 m	1/pk			
MANUFACTURE NAME	63 Community Drive Sanford, ME 04073 Phone: +1.866.736. Email: CustomerCar	GVS North America, Inc. 63 Community Drive Sanford, ME 04073 Phone: +1.866.736.1250 Email: CustomerCareGVSLS@gvs.com Website: www.gvs.com					
INTENDED USE/ APPLICATION	Typical Applications:  Protein and enzyme filtration and sterilization Biological fluid filtration and sterilization Pharmaceutical sterilization Environmental water studies Industrial water filtration IV filters IV drug delivery TPN delivery						
MATERIALS	Film: Polyethersulfone (PES), supported with polyester film						
	Regulatory Documentation Required:  □ Biocompatibility according ISO 10993-1  □ IMDS  □ Rohs, directive 2002/95/CE  □ BSE/TSE, directive 2003/32/CE  □ DEHP plasticizer and latex free  □ Aging  □ Reach 1907/2006/CE (hazardous substances regulation)  □ Dir. 67/548/CE and Reg. 1272/2008/CE (medical sector dangerous substances)  □ Conflict minerals  □ USP Class VI						



## PRODUCT SPECIFICATION

Product PN	3032875, 3032876, 3029505, 3018505, 1239465, 1235748, 3057106	Mod. 984
Description	0.03 µm Polyethersulfone Filtration Membrane Discs and Sheets	Rev. 01

TIOUUCI FIN 30320	75, 5052070, 5029505, 5	010303, 1239403, 1233740, 3037100		- WOU. 904
Description 0.03 μ	m Polyethersulfone Filt	ration Membrane Discs and Sheets		Rev. 01
PRODUCT CHARACTERISTICS	Low dr     Superior	ohilic ktractables ug and protein binding or burst strength ot consistency		
PRODUCT SHELF LIFE	5 years in normal stora Store the product at roo	age conditions: om temperate, in a dry place, free of dus	st and not exposed to dire	ect sunlight.
STERILIZATION	<ul> <li>This is a non-sterile pro</li> <li>⊠EtO</li> <li>⊠Gamma (25 kGy)6.5</li> <li>⊠Autoclave (121°C)</li> <li>*user has the responsibility to</li> </ul>	oduct, compatible with these sterilization  o validate their process	techniques*:	
PACKAGING AND LABELING	<ul><li>Quanti</li><li>Product</li><li>Manufa</li><li>Lot nur</li></ul>	ot item code and description acturing date	Exterior package Clear/white plastic box Round container White mailer  Cardboard box	
	Example label below:			

#### Example label below:



#### CERTIFICATE OF COMPLIANCE

A conformity declaration will be sent to the customer for each shipment with lot number and date of manufacture. The Quality management system is in compliance with ISO 9001.



## PRODUCT SPECIFICATION

Product PN	3032875, 3032876, 3029505, 3018505, 1239465, 1235748, 3057106	Mod. 984
Description	0.03 µm Polyethersulfone Filtration Membrane Discs and Sheets	Rev. 01

PERFORMANCE REQUIREMENTS		Acceptance Requirement		Standard Test Method
	1	Thickness	110-150 μm	N/A
	2	Water Flow Rate	Min 3 mL/min/cm <sup>2</sup> @ 10 psi	ASTM F 316-03 (2011)
	3	Bubble Point	90-110 psi (6.2-7.5 bar)	ASTM F 316-03 (2011)

This material specification describes the properties of product above indicated.

This document contains general requirements, material description, drawing references, defect specification, and biological material requirements.

#### **REVISIONS AND APPROVALS:**

DATE	REV	REASON FOR CHANGE	ISSUED AND CONTROLLED BY: (name /function and signature)	APPROVED BY: (name /function and signature)
6/8/18	0	Initial Release	Gil Tavares/Process Development Engineer	Kevin Wrigley/Director of Quality  Serm M. W.

### **CUSTOMER APPROVAL:**

We accept this n	naterial specification as a part of the agreed	terms of delivery.
Company name		
Approved by:		
	(Name, Function)	(Signature)
Date(Company stamp)		

Please send back this document signed for approval. If we will not receive this specification signed, we consider the first order placed as implicit approval.